



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

The Amsco Century Small Sterilizers are designed for sterilization of materials used in healthcare facilities. They are available in two configurations:

Gravity – designed for sterilization of heat and moisture-stable materials and sterilization of items intended for immediate use. Gravity sterilizers are equipped with gravity and flash cycles.

Prevacuum – designed for fast, efficient sterilization of heat and moisture-stable materials in addition to the same sterilization capabilities as a gravity sterilizer. Prevacuum sterilizers are equipped with prevacuum, gravity, flash, express, leak test and daily air removal test cycles.

Primary Product Features:

Control System

- Touch-sensitive screen with easy to read character display area.
- Ink on paper impact printer.
- Help screens for programming and trouble shooting alarm conditions.
- Communication interface compatible with most PC peripheral devices.
- Automatic check of control program and cycle data maintains process integrity.
- Service reprogrammable flash ROM memory.

Vertical Sliding Door

- Hands free loading and unloading.
- Foot Pedal activated door opening and closing.
- Non lubricated, steam-activated door seal

Modular Vessel and Piping

- Increased dependability and reduced service times.
- Reduced piping components for increased reliability.
- Emergency manual exhaust valve.
- Electronic water conservation valve.



Specifications:

Interior Chamber Dimensions

V116 16 x 16 x 26 inches
V120 20 x 20 x 38 inches

Product Configurations:

| | |
|--------------|-----------------------------|
| 16 x 16 x 26 | Gravity |
| 16 x 16 x 26 | Gravity with Liquid Cycle |
| 20 x 20 x 38 | Gravity |
| 20 x 20 x 30 | Gravity with Liquid Cycle |
| 16 x 16 x 26 | Prevacuum |
| 16 x 16 x 26 | Prevacuum with Liquid Cycle |
| 20 x 20 x 38 | Prevacuum |
| 20 x 20 x 30 | Prevacuum with Liquid Cycle |

Single Door Mounting

Cabinet Enclosed/Freestanding
Recessed

Double Door Mounting

Recessed through One Wall
Recessed through Two Walls**

Special Considerations:

16 x 16 x 26 double door sterilizers are not available with electric steam.

**Available for 20 x 20 x 38 inch double door sterilizers only.