



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

## Product Overview

The Aestiva 5 offers a good combination of quality and performance at an affordable price. The 7900 ventilator offer a full range of operational parameters including optimal tidal volume compensation, integrated volume and oxygen monitoring. The machine also features Tek 5 Vaporizers.

The Aestiva 5 has been optimized for patient safety through patient centered engineering while also being ergonomically configured for the clinician's comfort.



## Specifications:

Dimensions:	2 Vaporizers	3 Vaporizers	Absorber Bag Arms	Adjustable	Non-Adjustable
Height:	53.4 inches	53.4 inches	Arm Length:	12 inches	10 inches
Width:	29.5 inches	36.6 inches	Bag Arm Height:	34.3 inches	36 inches
Depth:	32.7 inches	32.7 inches		41 inches	
Weight:	300 lbs	340 lbs			

	Height	Width	Depth	Weight Limit
<b>Work Surface</b>	34.5 inches	18.5 inches	12.4 inches	N/A
<b>Folding Side Shelf (Optional)</b>	34.5 inches	10.5 inches	12.4 inches	25 lbs
<b>Top Drawer (1 std)</b>	4.1 inches	15.2 inches	10.2 inches	N/A
<b>Lower Drawers (Optional)</b>	5.7 inches	15.2 inches	10.2 inches	N/A

Top Shelves (Optional)	Width	Depth	Weight Limit
2 Vaporizers	18.7, 26.6 or 34.4 in	16.1 inches	100 lbs
3 Vaporizers	34.4 or 26.6 inches	16.1 inches	100 lbs

**Power Input** 120 VAC, 60Hz, 220-240 VAC 50Hz, 230 VAC 50Hz, 100 VAC 50/60 Hz, 240 VAC 50 Hz.  
**Absorber** Rotation 850  
**Castors** Two Locking Front Castors