



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

Product Overview

The Datex Ohmeda Modulus SE Anesthesia Machine features a Selecta-Tek Back Bar that facilitates mounting 3 Tek-5 vaporizers. It is supplied standard with the model 7800 electronic ventilator. Each unit is supplied standard with a 3-drawer module and overhead shelving. Like all Ohmeda machines the supply pressure gauges, flow meter controls and ventilator functions have all been located within easy reach and in clear view of the care giver.

Each machine is supplied standard with one vaporizer (Isoflurane or Halothane). Gas supply hoses should be ordered separately.



**DATEX OHMEDA MODULUS SE
ANESTHESIA MACHINE**

Specifications:

Monitoring Includes: Inspired Volume, Pressure and Oxygen Concentration

Key Features: Ascending bellows ventilator (Model 7800). Manual/Spontaneous breathing, mechanical flow meters (rotameter/flow tube) and vaporizers. Traditional size and configuration of absorber head and higher volume breathing circuit. Loose absorbent granules or canisters. Removable vaporizers mounted on a Selecta Tek Back Bar. Can be used with non-rebreathing circuits. Opened or closed scavenger interface.

Machine Weight: 306 lbs