



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

## Product Overview

The most accurate delivery of anesthetic agents over a flow range of 300 mls to 15L, which provides excellent reliability.

### Options:

- Forane 19.1 vaporizer
- Halothane 19.1 vaporizer
- Ethrane 19.1 vaporizer



## DRAGER 19.1 VAPORIZERS

### Specifications:

**Weight:**.....6 kg  
**Height:**.....205 mm  
**Diameter:**.....100 mm  
**Temperature range:**.....10 - 40 degrees Celsius  
**Flow Range:**.....0.3 - 15l/min

### Pressure Compensation Range:

Pressure difference at pressure change:.....20 mbar  
Maximum Pressure:.....200 mbar  
Maximum Inclination:.....45 degrees