



3915 152ND ST NE | MARYSVILLE, WA 98271 | 1.800.827.3747

## Microscope Overview

Notable Features of the OPMI-6 are:

- Motorized Focus and Zoom
- New Fiber-Optic Light Source and Cable
- Foot pedal for Motorized Features
- Various Optics Available – OPMI 6, 6CFC, 6S

Spare Bulbs, Sterilizable Knob Covers and Dust Cover provided with Every Microscope System.

2 Year Warranty Including Parts and Labor.



## ZEISS OPMI - 6 SERIES

### Specifications:

**Binocular:**.....Inclined Inclinable Upgrade Available

**Objective Lens:**.....Your Choice: 150mm-400mm

**Magnification :**..... 1:4 Ratio; Motorized

**Focused:**.....Motorized, 50mm Range

**Eyepieces:**.....12.5X HEP

**Pupillary Distance:**.....50mm-70mm

**Filter:**.....Standard UV

**Light Source:**.....Variable 15V, 150W

**Illuminated Field:**.....Up to 48mm

**Power Source:**..... 110V-240V; 50-60 Hz

**Mounting System:**..... Floorstand w/Locking Castors

**Horizontal Reach:**..... Arm 24 inches

**End of Optics:**..... 29 inches

**Vertical Range (Macro Focus):**.....65 inches

**Arm:**.....Pole Style, Sliding

**Shipping Weight:**.....200 lbs

**Warranty Term:**.....2 Years Parts and Labor